


## PMU DEVICES & ACCESSORIES HYGIENE CARTRIDGE DESCRIPTION


For the development of the various hygiene cartridges we use only the latest technology, high-quality materials and our many years of experience. We know that every PMU stylist develops their own style over time. That is why a wide range of different

hygiene cartridges is an advantage. The basis of the hygiene cartridges is a medically approved plastic, the needles are made of stainless steel. To ensure perfect quality, the production of our cartridges is kept under constant and independent control.




**1 – NANO NT**  
Fine lip contours, fine hairs; needle: stiff, not flexible, no vibration, good for beginners.  
Ø 0.25 mm

**G C**




**1 – NANO (1)**  
Fine hairs, fine lip contours, capillary drawing on areola, fine shading techniques; needle: flexible, very sharp.  
Ø 0.20 mm

**G C S T**




**1 – NANO (2)**  
PMU removal, fine hairs, fine lip contouring, capillary drawing on areola, fine shading techniques (pixels); Needle: medium flexible.  
Ø 0.25 mm

**C S**




**1 – MICRO**  
Exact lip contouring, hair and shading shading techniques, head pigmentation; Needle: stiff.  
Ø 0.4 mm

**G C S T**




**1 – LINER**  
Intensive, wide eyelid lines, eyebrow shading techniques, intensive lip contouring, scalp pigmentation.  
Ø 0.4 mm

**G C S T**




**3 – MICRO**  
PMU removal, eyeliner, lash line, eyebrow shading, lip contouring, lip shading, areola, scars.  
Ø 0.18 mm

**G C T**




**3 – OUTLINE**  
PMU removal, eyeliner, lash line, eyebrow shading (Pixel/Powder technique), lip contouring, lip shading.  
Ø 0.25 mm

**G C S T**




**3 – LINER**  
Eyebrow shading techniques, Lip shading, wide lip contour.  
Ø 0.3 mm

**G C S T**




**3 – POWER**  
Compact, intensive shading techniques: Lips, areola, eyeliner, lash line, lip contour.  
Ø 0.3 mm

**C T**




**3 – SLOPE**  
Fine hairs, techniques for shading the eyebrows, eyeliner with shading, lash line, lip contour, lip shading, areola.  
Ø 0.4 mm

**C S T**




**3 – NANO SLOPE**  
Digital blading for fine hairs, light regular shading techniques for: eyebrows, eye shadow, full lips, areola.  
Ø 0.25 mm

**G C**



**3 – NANO SPREAD**  
Pixel/Powder shading techniques for: eyebrows, eyeliner, lips, areola; dermaneedling of small scars; needles: flexible.  
Ø 0.25 mm


**G**



**4 – FLAT**  
Dermaneedling of stretch marks, hair thickening; shading techniques: light shadows, scalp pigmentation, lips.  
Ø 0.35 mm


**C S**

- G** Genius hygiene cartridges
- C** Classic hygiene cartridges
- S** Symphony hygiene cartridges
- T** TOP hygiene cartridges




**5 – ROUND**  
Shading techniques: Lips, eyebrows, areola, PMU removal.  
Ø 0.3 mm

**C**




**5 – SHADER**  
PMU removal; Pixel/Powder shading technique, eyebrows, lips.  
Ø 0.3 mm

**G C S T**




**5 – POWER**  
Shading technique Pixel/Powder, areola, lips, dermaneedling.  
Ø 0.3 mm

**C**




**5 – MAGNUM**  
Shading techniques: Areola, lips, dermaneedling of scars, hair thickening, large areas.  
Ø 0.3 mm

**C S**




**5 – SLOPE**  
Shading techniques: areola, lips, eyebrows.  
Ø 0.4 mm

**C S**



**5 – NANO SLOPE**  
Digital blading for fine hairs, light regular shading techniques for eyebrows, eyeliner with shading, lip shading, areola.  
Ø 0.25 mm


**G C**



**5 – NANO SPREAD**  
Shading techniques Pixel/Powder for: eyebrows, eyeliner, lips, areola, dermaneedling; scars of all kinds, large areas; needles: flexible.  
Ø 0.25 mm


**G**

## PMU DEVICES & ACCESSORIES HYGIENE CARTRIDGE DESCRIPTION




**7 – POWER**  
Full shading techniques: Areola, lips, dermaneedling.  
Ø 0.3 mm

**C**




**7 – ROUND**  
PMU removal, full shading techniques: areola, lips, dermaneedling. Scars of all kinds, large areas.  
Ø 0.3 mm

**C S**




**7 – NANO ROUND**  
Shading techniques: areola, lips; dermaneedling, hair thickening; needles: flexible.  
Ø 0.25 mm

**G**




**9 – MAGNUM**  
Shading techniques: areola, lips; dermaneedling, hair thickening, scars of all kinds on large areas.  
Ø 0.3 mm

**C S**



**DERMA (6)**  
Dermaneedling for scars, wrinkles, hair thickening; needles are fixed to the movable module head.  
Ø 0.3 mm

**C**



**DERMA (18)**  
Dermaneedling of wrinkles, pigment spots, décolleté, neck; needles are fixed to the movable module head.  
Ø 0.3 mm




**C S T**




# PMU DEVICES & ACCESSORIES









## HYGIENE CARTRIDGE OVERVIEW

- G** Genius hygiene cartridges
- C** Classic hygiene cartridges
- S** Symphony hygiene cartridges
- T** TOP hygiene cartridges

-  Lips
-  Eyeliner
-  Eyebrow

-  Correction
-  Areola
-  Camouflage, Face

-  Scalp,
-  Hair thickening
-  Dermeneedling

	GROUND	CARTRIDGE OVERVIEW				TREATMENT OPTIONS							
	S/M/L	<b>G</b>	<b>C</b>	<b>S</b>	<b>T</b>								
1 – NANO NT	S	•	•			•		•					
1 – NANO (1)	S	•	•	•	•	•		•		•			
1 – NANO (2)	S		•	•		•		•	•	•			
1 – MICRO	L	•	•	•	•	•	•	•				•	
1 – LINER	M	•	•	•	•	•	•	•				•	
3 – MICRO	S	•	•		•	•	•	•	•	•	•		
3 – OUTLINE	M	•	•	•	•	•	•	•	•		•		
3 – LINER	M	•	•	•	•	•	•	•					
3 – POWER	L		•		•	•	•			•			
3 – SLOPE	M		•	•	•	•	•	•		•			
3 – NANO SLOPE	M	•	•			•	•	•		•			
3 – NANO SPREAD	S	•				•	•	•		•			•
4 – FLAT	M		•	•		•						•	•
5 – ROUND	M		•			•		•	•	•			
5 – SHADER	M	•	•	•	•	•		•	•				
5 – POWER	L		•			•					•		•
5 – MAGNUM	M		•	•		•				•	•		•
5 – SLOPE	M		•	•		•		•		•			
5 – NANO SLOPE	S	•	•			•	•	•		•			
5 – NANO SPREAD	S	•				•		•		•	•		•
7 – POWER	L		•			•				•			•
7 – ROUND	M		•	•				•		•	•		•
7 – NANO ROUND	S	•				•				•	•	•	•
9 – MAGNUM	M		•	•						•	•	•	•
DERMA (6)	M		•								•	•	•
DERMA (18)	M		•	•	•								•

- S** Short milled/ground : fast pigment flow, shading: 80 rpm, contouring: 90 rpm, hand speed: fast to medium
- M** Medium milled/ground: gentle pigment flow, shading: 70 rpm, contouring: 90 rpm, hand speed: medium
- L** Long milled/ground: slow pigment flow, max. 90 rpm, hand speed: smooth, slow