














## SWISS COLOR® PMU DEVICES & ACCESSORIES HYGIENE CARTRIDGE DESCRIPTION













For the development of the various hygiene cartridges we use only the latest technology, high quality materials and our many years of experience. We know that every PMU stylist develops her/his own style over time. That is why a wide range of different

hygiene cartridges is an advantage. The basis of the Hygiene cartridges is a medically approved plastic, and the needles are made of stainless steel. To ensure the perfect quality, the production of our cartridges is kept under constant and independent control.

 <p><b>1 – NANO NT</b> Fine lip outline, fine hairs; needle: stiff, not flexible, no vibration, good for beginners. Ø 0.25 mm</p> <p><b>G C</b></p>	 <p><b>3 – LINER</b> Eyebrow shading techniques, lip shading, thicker lip contours. Ø 0.3 mm</p> <p><b>G C S T</b></p>
 <p><b>1 – NANO (1)</b> Fine hairs, fine lip contours, capillary drawing on areola, fine shading techniques; needle: flexible, very sharp. Ø 0.20 mm</p> <p><b>G C S T</b></p>	 <p><b>3 – POWER</b> Compact, intensive shading techniques: lips, areola, eyeliner, lash line, lip contour. Ø 0.3 mm</p> <p><b>C T</b></p>
 <p><b>1 – NANO (2)</b> PMU removal, fine hairs, fine lip contours, capillary drawing on areola, fine shading techniques (pixels); needle: medium flexible. Ø 0.25 mm</p> <p><b>C S</b></p>	 <p><b>3 – SLOPE</b> Fine hairs, techniques for shading the eyebrows, eyeliner with shading, lash line, lip contour, lip shading, areola. Ø 0.4 mm</p> <p><b>C S T</b></p>
 <p><b>1 – MICRO</b> Exact lip contouring, hairs and shading techniques, scalp pigmentation; needle: stiff. Ø 0.4 mm</p> <p><b>G C S T</b></p>	 <p><b>3 – NANO SLOPE</b> Digital blading for fine hairs, light regular shading techniques for: eyebrows, eye shadow, full lips, areola. Ø 0.25 mm</p> <p><b>G C</b></p>
 <p><b>1 – LINER</b> Intense, wide eyelid lines, eyebrow shading techniques, intense lip contouring, scalp pigmentation. Ø 0.4 mm</p> <p><b>G C S T</b></p>	 <p><b>3 – NANO SPREAD</b> Pixel/Powder shading techniques for: eyebrows, eyeliner, lips, areola; dermaneedling of small scars; needles: flexible. Ø 0.25 mm</p> <p><b>G</b></p>
 <p><b>3 – MICRO</b> PMU removal, eyeliner, lash line, eyebrow shading, lip contouring, lip shading, areola, scars. Ø 0.18 mm</p> <p><b>G C T</b></p>	 <p><b>4 – FLAT</b> Dermaneedling of stretch marks, hair thickening; shading techniques: light shadows, scalp pigmentation, lips. Ø 0.35 mm</p> <p><b>C S T</b></p>
 <p><b>3 – OUTLINE</b> PMU removal, eyeliner, lash line, eyebrow shading (pixel/powder technique), lip liner, lip shading. Ø 0.25 mm</p> <p><b>G C S T</b></p>	

## SWISS COLOR® PMU DEVICES & ACCESSORIES HYGIENE CARTRIDGE DESCRIPTION

**G** SC Genius PRO 2.5 Hygiene cartridges  
**C** Classic Hygiene cartridges  
**S** Symphony Hygiene cartridges  
**T** TOP Hygiene cartridges

 <p><b>5 – ROUND</b> Shading techniques: lips, eyebrows, areola, PMU removal. Ø 0.3 mm</p> <p><b>C T</b></p>	 <p><b>7 – POWER</b> Full shading techniques: areola, lips, dermaneedling. Ø 0.3 mm</p> <p><b>C</b></p>
 <p><b>5 – SHADER</b> PMU removal; Pixel/Powder shading technique, eyebrows, lips. Ø 0.3 mm</p> <p><b>G C S T</b></p>	 <p><b>7 – ROUND</b> PMU removal, full shading techniques: areola, lips, dermaneedling. Scars of all types, large areas. Ø 0.3 mm</p> <p><b>C S T</b></p>
 <p><b>5 – POWER</b> Shading technique Pixel/Powder, areola, lips, dermaneedling. Ø 0.3 mm</p> <p><b>C</b></p>	 <p><b>7 – NANO ROUND</b> Shading techniques: areola, lips; dermaneedling, hair thickening; needles: flexible. Ø 0.25 mm</p> <p><b>G</b></p>
 <p><b>5 – MAGNUM</b> Shading techniques: areola, lips, dermaneedling of scars, hair thickening, large areas. Ø 0.3 mm</p> <p><b>C S</b></p>	 <p><b>9 – MAGNUM</b> Shading techniques: areola, lips; dermaneedling, hair thickening, scars of all kinds on large areas. Ø 0.3 mm</p> <p><b>C S T</b></p>
 <p><b>5 – SLOPE</b> Shading techniques: areola, lips, eyebrows. Ø 0.4 mm</p> <p><b>C S</b></p>	 <p><b>DERMA (6)</b> Dermaneedling for scars, wrinkles, hair thickening; needles are fixed on the movable cartridge head. Ø 0.3 mm</p> <p><b>C</b></p>
 <p><b>5 – NANO SLOPE</b> Digital blading for fine hairs, light regular shading techniques for eyebrows, eyeliner with shading, lip shading, areola. Ø 0.25 mm</p> <p><b>G C</b></p>	 <p><b>DERMA (18)</b> Dermaneedling for wrinkles, pigment irritations, décollete, throat; needles are fixed on the movable cartridge head. Ø 0.3 mm</p> <p><b>C S T</b></p>
 <p><b>5 – NANO SPREAD</b> Shading technique Pixel/Powder for: eyebrows, eyeliner, lips, areola, dermaneedling; scars of all kinds, large areas; needles: flexible. Ø 0.25 mm</p> <p><b>G</b></p>	



SWISS COLOR® PMU DEVICES & ACCESSORIES  
HYGIENE CARTRIDGE OVERVIEW

	GROUND	CARTRIDGE OVERVIEW				TREATMENT OPTIONS							
	S/M/L	G	C	S	T	Lips	Eye	Brow	Hand	Areola	Camouflage	Scalp	Derma
1 – NANO NT	S	●	●			●		●					
1 – NANO (1)	S	●	●	●	●	●		●		●			
1 – NANO (2)	S		●	●		●		●	●	●			
1 – MICRO	L	●	●	●	●	●	●	●				●	
1 – LINER	M	●	●	●	●	●	●	●				●	
3 – MICRO	S	●	●		●	●	●	●	●	●	●		
3 – OUTLINE	M	●	●	●	●	●	●	●	●		●		
3 – LINER	M	●	●	●	●	●	●	●					
3 – POWER	L		●		●	●	●			●			
3 – SLOPE	M		●	●	●	●	●	●		●			
3 – NANO SLOPE	M	●	●			●	●	●		●			
3 – NANO SPREAD	S	●				●	●	●		●			●
4 – FLAT	M		●	●	●	●						●	●
5 – ROUND	M		●		●	●		●	●	●			
5 – SHADER	M	●	●	●	●	●		●	●				
5 – POWER	L		●			●					●		●
5 – MAGNUM	M		●	●		●				●	●		●
5 – SLOPE	M		●	●		●		●		●			
5 – NANO SLOPE	S	●	●			●	●	●		●			
5 – NANO SPREAD	S	●				●		●		●	●		●
7 – POWER	L		●			●				●			●
7 – ROUND	M		●	●	●	●		●		●	●		●
7 – NANO ROUND	S	●				●				●	●	●	●
9 – MAGNUM	M		●	●	●	●				●	●	●	●
DERMA (6)	M		●								●	●	●
DERMA (18)	M		●	●	●								●

**S** Short milled/ground : fast pigment flow, shading: 80 rpm, contouring: 90 rpm, hand speed: fast to medium  
**M** Medium milled/ground: gentle pigment flow, shading: 70 rpm, contouring: 90 rpm, hand speed: medium  
**L** Long milled/ground: slow pigment flow, max. 90 rpm, hand speed: smooth, slow

**G** SC Genius PRO 2.5 Hygiene cartridges  
**C** Classic Hygiene cartridges  
**S** Symphony Hygiene cartridges  
**T** TOP Hygiene cartridges

Lips  
 Eyeliner  
 Eyebrow

Correction  
 Areola  
 Camouflage, Face

Scalp,  
 Hair thickening  
 Dermaneedling