


PMU DEVICES & ACCESSORIES HYGIENE CARTRIDGE DESCRIPTION


For the development of the various hygiene cartridges we use only the latest technology, high-quality materials and our many years of experience. We know that every PMU stylist develops their own style over time. That is why a wide range of different

hygiene cartridges is an advantage. The basis of the hygiene cartridges is a medically approved plastic, the needles are made of stainless steel. To ensure perfect quality, the production of our cartridges is kept under constant and independent control.




1 – NANO NT
Fine lip contours, fine hairs; needle: stiff, not flexible, no vibration, good for beginners.
Ø 0.25 mm

G C




1 – NANO (1)
Fine hairs, fine lip contours, capillary drawing on areola, fine shading techniques; needle: flexible, very sharp.
Ø 0.20 mm

G C S T




1 – NANO (2)
PMU removal, fine hairs, fine lip contouring, capillary drawing on areola, fine shading techniques (pixels); needle: medium flexible.
Ø 0.25 mm

C S




1 – MICRO
Exact lip contouring, hair and shading techniques, head pigmentation; needle: stiff.
Ø 0.4 mm

G C S T




1 – LINER
Intensive, wide eyelid lines, eyebrow shading techniques, intensive lip contouring, scalp pigmentation.
Ø 0.4 mm

G C S T




1 – LINER CAPILLARY
Fine hairs and contours, pixel shading; eyebrows, eyeliner, lips; more pigment absorption and improved color flow for longer working time.
Ø 0.4 mm

G




3 – MICRO
PMU removal, eyeliner, lash line, eyebrow shading, lip contouring, lip shading, areola, scars.
Ø 0.18 mm

G C T




3 – OUTLINE
PMU removal, eyeliner, lash line, eyebrow shading (Pixel/Powder technique), lip contouring, lip shading.
Ø 0.25 mm

G C S T




3 – LINER
Eyebrow shading techniques, lip shading, wide lip contour.
Ø 0.3 mm

G C S T




3 – POWER
Compact, intensive shading techniques: lips, areola, eyeliner, lash line, lip contour.
Ø 0.3 mm

C T




3 – SLOPE
Fine hairs, techniques for shading the eyebrows, eyeliner with shading, lash line, lip contour, lip shading, areola.
Ø 0.4 mm

C S T




3 – NANO SLOPE
Digital blading for fine hairs, light regular shading techniques for: eyebrows, eye shadow, full lips, areola.
Ø 0.25 mm

G C



3 – NANO SPREAD
Pixel/Powder shading techniques for: eyebrows, eyeliner, lips, areola; dermanedling of small scars; needles: flexible.
Ø 0.25 mm


G



4 – FLAT
Dermaneedling of stretch marks, hair thickening; shading techniques: light shadows, scalp pigmentation, lips.
Ø 0.35 mm


C S

- G** Genius hygiene cartridges
- C** Classic hygiene cartridges
- S** Symphony hygiene cartridges
- T** TOP hygiene cartridges




5 – ROUND
Shading techniques: lips, eyebrows, areola, PMU removal.
Ø 0.3 mm

C




5 – SHADER
PMU removal; Pixel/Powder shading technique, eyebrows, lips.
Ø 0.3 mm

G C S T




5 – POWER
Shading technique Pixel/Powder, areola, lips, dermaneedling.
Ø 0.3 mm

C




5 – MAGNUM
Shading techniques: areola, lips, dermaneedling of scars, hair thickening, large areas.
Ø 0.3 mm

C S




5 – SLOPE
Shading techniques: areola, lips, eyebrows.
Ø 0.4 mm

C




5 – NANO SLOPE
Digital blading for fine hairs, light regular shading techniques for eyebrows, eyeliner with shading, lip shading, areola.
Ø 0.25 mm

G C




5 – NANO SPREAD
Shading techniques Pixel/Powder for: eyebrows, eyeliner, lips, areola, dermaneedling; scars of all kinds, large areas; needles: flexible.
Ø 0.25 mm

G




7 – POWER
Full shading techniques: areola, lips, dermaneedling.
Ø 0.3 mm

C




7 – ROUND
PMU removal, full shading techniques: areola, lips, dermaneedling. Scars of all kinds, large areas.
Ø 0.3 mm

C S




7 – NANO ROUND
Shading techniques: areola, lips; dermaneedling, hair thickening; needles: flexible.
Ø 0.25 mm

G




9 – MAGNUM
Shading techniques: areola, lips; dermaneedling, hair thickening, scars of all kinds on large areas.
Ø 0.3 mm

C S



DERMA (6)
Dermaneedling for scars, wrinkles, hair thickening; needles are fixed to the movable module head.
Ø 0.3 mm

C



DERMA (18)
Dermaneedling of wrinkles, pigment spots, décolleté, neck; needles are fixed to the movable module head.
Ø 0.3 mm

C S T




PMU DEVICES & ACCESSORIES HYGIENE CARTRIDGE DESCRIPTION




PMU DEVICES & ACCESSORIES









HYGIENE CARTRIDGE OVERVIEW

- G** Genius hygiene cartridges
- C** Classic hygiene cartridges
- S** Symphony hygiene cartridges
- T** TOP hygiene cartridges

-  Lips
-  Eyeliner
-  Eyebrow

-  Correction
-  Areola
-  Camouflage, Face

-  Scalp,
-  Hair thickening
-  Dermeneedling

	GROUND	CARTRIDGE OVERVIEW				TREATMENT OPTIONS							
	S/M/L	G	C	S	T								
1 – NANO NT	S	•	•			•		•					
1 – NANO (1)	S	•	•	•	•	•		•		•			
1 – NANO (2)	S		•	•		•		•	•	•			
1 – MICRO	L	•	•	•	•	•	•	•				•	
1 – LINER	M	•	•	•	•	•	•	•				•	
1 – LINER CAPILLARY	S	•				•	•	•					
3 – MICRO	S	•	•		•	•	•	•	•	•	•		
3 – OUTLINE	M	•	•	•	•	•	•	•	•		•		
3 – LINER	M	•	•	•	•	•	•	•					
3 – POWER	L		•		•	•	•			•			
3 – SLOPE	M		•	•	•	•	•	•		•			
3 – NANO SLOPE	M	•	•			•	•	•		•			
3 – NANO SPREAD	S	•				•	•	•		•			•
4 – FLAT	M		•	•		•						•	•
5 – ROUND	M		•			•		•	•	•			
5 – SHADER	M	•	•	•	•	•		•	•				
5 – POWER	L		•			•					•		•
5 – MAGNUM	M		•	•		•				•	•		•
5 – SLOPE	M		•			•		•		•			
5 – NANO SLOPE	S	•	•			•	•	•		•			
5 – NANO SPREAD	S	•				•		•		•	•		•
7 – POWER	L		•			•				•			•
7 – ROUND	M		•	•				•		•	•		•
7 – NANO ROUND	S	•				•				•	•	•	•
9 – MAGNUM	M		•	•						•	•	•	•
DERMA (6)	M		•								•	•	•
DERMA (18)	M		•	•	•								•

- S** Short milled/ground : fast pigment flow, shading: 80 rpm, contouring: 90 rpm, hand speed: fast to medium
- M** Medium milled/ground: gentle pigment flow, shading: 70 rpm, contouring: 90 rpm, hand speed: medium
- L** Long milled/ground: slow pigment flow, max. 90 rpm, hand speed: smooth, slow